

RX3 Vortex Mixer

Vortex Mixer for Basic Requirements

- Fixed stirring speed, 3000 rpm
- Touch mode
- Enhanced stability on the bench
- New materials, higher performance



The RX3 is a basic vortex mixer, operating in **touch mode**.

High-speed and excellent vortexing is ensured by the fixed speed at 3000 rpm. The RX3 remains **extremely stable** thanks to special feet.

Get full speed by pressing down the rubber cup and mix your sample in an excellent way.

The new technopolymer structure provides an **optimum chemical resistance**.

The RX3 is ideal for many applications, mixing samples in different containers.

VELP Scientifica offers many years of high performance in laboratories all over the world, with resistant and reliable solutions.

Live an amazing vortex experience!

Features and Benefits

Optimized Performance, Quiet Operation

Fixed speed (3000 rpm).

Product Handling

Ergonomic and highly innovative design for excellent stability on the bench.

It is usable on many surfaces thanks to the 4 anti-sliding feet.

Technical Data	Description
Construction material:	technopolymer structure
Support system:	4 anti-sliding feet for strong fixing and high stability
Protection rating CEI EN 60529:	IP 42
Electronic speed regulation:	constant, 3000 rpm
Power Supply:	100-240V/50-60Hz
Power:	15 W
Weight:	2,7 Kg (5.9 lb)
Dimensions (W x H x D) :	150x130x165 mm (5.9x5.1x6.5 in)
Ordering Information Code No	Description
F202A0171	RX3 Vortex Mixer 100-240V/50-60Hz
A00000016	Small rubber supporting plate Ø 50 mm
10003083	US plug
10003084	UK plug
10003085	Australian plug



A00000016

Your authorized agent:

We reserve the right to make technical alternations
We do not assume liability for errors in printing, typing or transmission



VELP Scientifica srl
via Stazione 16
20040 Usmate (Milano) Italy
Tel +39 039 628811
Fax +39 039 6288120
inse@velp.it
www.velp.com